



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

cappel

Catalog Number: 10706
Anti-Synaptophysin Polyclonal

SPECIES: Rabbit

Ig CLASS: IgG

REAGENT: Serum

IMMUNOGEN: Synthetic human synaptophysin

SPECIFICITY: Reactive in human and rat

POSITIVE CONTROL: Pituitary gland

APPLICATIONS:

- Frozen sections; fixation: acetone, 5 minutes, -20°C.
- Paraffin sections (no pretreatment); fixation: formalin.
- Cell Preparations

DILUTION: 1:50-1:100 (Indirect method, 30 minute incubation). Dilute immediately before use with PBS. Optimal dilution should be tested by serial dilution.

REFERENCE:

- Luts, L., et al. **Anat. Embryol.** 195: 515-524 (1997).

Note: This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.

[If a titer/working dilution is not given above, please click here to see a general dilution chart for working with antibodies. Please note that the general dilution chart should only be used as a guideline. Each lab should determine their own optimal working dilution.](#)

[Will this antibody work with your application? Please click here to see a general chart of antibody applications. Please note that any information given above is primary application data. The general applications charts should only be used as a reference.](#)