



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

cappel

Catalog Number: 55838

Purified Human Gamma Globulin

DESCRIPTION: Product is the lyophilized powder of human gamma globulin and buffer salts.

Human gamma globulin is purified from pooled human serum using multi-step procedures which may include salt fractionation, ion-exchange chromatography and immunoabsorption. The product is dialyzed into 0.01M sodium phosphate, with 0.07M sodium chloride, pH 7.3, filtered through a 0.22 µm filter, vialled and lyophilized.

USES: This product is suitable for use as a carrier protein in secondary antibody systems, e.g., radioimmunoassay. Purified proteins are recommended as ELISA coating antigens and as controls.

See **Related Product Ordering Information**.

SPECIFICATIONS: The product is tested for purity at a minimum of 20 mg/ml using immunoelectrophoresis. The product is primarily human IgG, IgA, IgM.

The total protein is determined using the Biuret procedure, with bovine albumin as the standard. Each vial is over filled by approximately 10 per cent so that a minimum of 100 mg of human gamma globulin is obtained upon reconstitution.

RECONSTITUTION: Reconstitute product with 6 ml of deionized or distilled water. Gentle swirling may be used to speed rehydration. Avoid vigorous shaking of the constituted material.

WARNING AND CAUTION:

All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

RELATED PRODUCT ORDERING INFORMATION:

Catalog #DescriptionSize

55127 GOAT ANTI-HUMAN GAMMA GLOBULIN5ML
55100 GOAT ANTI-HUMAN IMMUNOGLOBULINS (G,A,M)2ML
55087 GOAT ANTI-HUMAN IGG (WHOLE MOLECULE) 2ML
55082 GOAT ANTI-HUMAN IGA (Å CHAIN)2ML
55097 GOAT ANTI-HUMAN IGM (5 FC µ)2ML
55020 IGG GOAT ANTI-HUMAN IMMUNOGLOBULINS (G,A,M) 5ML
55007 IGG GOAT ANTI-HUMAN IGG (WHOLE MOLECULE)5ML
55003 IGG GOAT ANTI-HUMAN IGA (Å CHAIN)5ML
55017 IGG GOAT ANTI-HUMAN IGM (5FC µ)5ML
55908 PURIFIED HUMAN IGG50MG
55906 PURIFIED HUMAN SERUM IGA 5MG
55916 PURIFIED HUMAN IGM5MG

Note: This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.