



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

cappel

Catalog Number: 57005

Fluorescein-Conjugated Rabbit Affinity Purified Antibody To Sheep IgG (Whole Molecule)

DESCRIPTION: Product is fluorescein-5-isothiocyanate (FITC "Isomer I")-conjugated rabbit affinity purified antibody to sheep IgG (whole molecule) and buffer salts.

Antibody and FITC are conjugated under defined conditions to obtain the specified molar ratio (3-6 mol/mol). Conjugated protein is purified by gel filtration. The product is dialyzed into 0.02M sodium phosphate, 0.14M sodium chloride, pH 7.3, with 10% glycerol (v/v), adjusted to standard titer, filtered through a 0.22 µm filter, and vialled.

USES: This product is suitable for use as a reagent in immunofluorescence assays (IFA), cell staining (fluorescent microscopy and cell sorting), tissue staining, and blot immunostaining. Use of these affinity purified antibodies is recommended to avoid non-specific binding from the inherent antibodies of host animals. (Note: F(ab')₂ fragments are recommended for staining of cells or tissues which contain Fc receptors. See **Related Product Ordering Information** for these products.)

SPECIFICATIONS: The total protein is measured by absorbance (A_{280nm}); F/P ratio is calculated using absorbance at 491.5nm. Antibody titer of the affinity purified antibody is standardized by equivalence-point precipitation. The affinity purified antibody is tested for purity at 10 mg/ml using immunoelectrophoresis. The product is rabbit IgG; no trace of albumin is detected. This product shows reactivity to sheep IgG; cross-reactivity to other species may exist. Because the product is directed against whole IgG, some antibodies are expected to bind light chain sites common to all immunoglobulins. Antibody activity to non-immunoglobulin serum proteins is not present.

RELATED PRODUCT ORDERING INFORMATION:

Catalog #DescriptionSize

55794 RABBIT ANTI-SHEEP IGG (WHOLE MOLECULE)2ML
55802 IGG RABBIT ANTI-SHEEP IGG (WHOLE MOLECULE) 5ML
55807 F(AB')₂ RABBIT ANTI-SHEEP IGG (WHOLE MOLECULE)2ML
55809 FITC RABBIT ANTI-SHEEP IGG (WHOLE MOLECULE) 2ML
55811 FITC F(AB')₂ RABBIT ANTI-SHEEP IGG (WM)2ML
55812 RHO RABBIT ANTI-SHEEP IGG (WHOLE MOLECULE)2ML
55814 HRP RABBIT ANTI-SHEEP IGG (WHOLE MOLECULE)2ML
56991AKP AFF.PUR. RABBIT ANTI-SHEEP IGG (WM)1ML
56992HRP AFF.PUR. RABBIT ANTI-SHEEP IGG (WM)2ML
55954 PURIFIED SHEEP IGG50MG
55990 NORMAL RABBIT SERUM2ML

Note: This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.

[If a titer/working dilution is not given above, please click here to see a general dilution chart for working with antibodies. Please note that the general dilution chart should only be used as a guideline. Each lab should determine their own optimal working dilution.](#)

[Will this antibody work with your application? Please click here to see a general chart of antibody applications. Please note that any information given above is primary application data. The general applications charts should only be used as a reference.](#)